

Extensive
Cloud based
Report

Digital Health
enabled

EMR/EHR
synchronized



Touchscreen
GUI Console


Completely
Stand-Alone

No Laptop/
PC needed


Non-Invasive Baro-Reflex Assessment – Combined System



 Four Parasympathetic Tests Based on Lead II ECG via a Stable 5-patch System

 Heart Rate Variability including exhaustive Time and Frequency Domain Analysis

 Two Sympathetic Tests based on Cuff BP

 All as per Ewings Tests recommended by the American Diabetes Association (ADA) as Gold Standard!



Supported by: